



Neuberg
DIAGNOSTICS

• India • UAE • South Africa • USA

Neu INSIGHTS



Primary Immunodeficiency Disorder

Serial number : 015 Edition : 1. 2022

What are Primary Immune Deficiencies (PIDs)?

- ▶ PIDs are a group of inherited disorders that affect the development and/or function of the immune system. As per an estimate, 1 in 1200 individuals are affected with a PID. This incidence however is increased in populations with a higher degree of consanguinity. These diseases usually result in recurrent or persistent severe infections, autoimmunity, auto-inflammation or malignancies. Though not uncommon, they are grossly underdiagnosed due to lack of awareness about these diseases.
- ▶ While most children with persistent infections may have a normal underlying immune system, it is important to recognize a child with an underlying PID so that further investigations can be performed. Prompt and early diagnosis not only helps initiate appropriate treatment, but also helps further genetic counselling for the family.




Primary Immunodeficiency
diseases affect more than

Million
people worldwide

With no gender, age or geographical
boundaries.

When to suspect PID?

Careful clinical evaluation is crucial for identification of PID. Ten warning signs have been suggested by the European society for Immunodeficiencies (ESID):




Patients showing these signs and symptoms must be subjected to further evaluation to confirm/exclude a PID.

Common clinical manifestations of the major groups of PIDs include

- ▶ **Combined T and B deficiency** – Opportunistic infections (Systemic viral infections, fungal and bacterial infections).
- ▶ **Antibody deficiency** – Upper and lower respiratory tract, GI tract, skin infections, sepsis and meningitis. Common organisms- encapsulated bacteria (Pneumococcus, H influenza)
- ▶ **Phagocytic defects** – suppurative infections – lung abscess, liver abscess, lymphadenitis, osteomyelitis. Organisms – S aureus, Klebsiella, Aspergillus, Candida, Burkholderia
- ▶ **Complement deficiency** – Infections with encapsulated bacteria eg: Meningococcal meningitis.

Clinical approach

Careful clinical evaluation is crucial for identification of PID. Ten warning signs have been suggested by the European society for Immunodeficiencies (ESID):




Screening tests for PIDs

Genetic tests are the standard diagnostic methods for most PIDs, but these are time consuming, expensive and not easily available. More easily available screening tests can be used to narrow down the possible diagnosis and then carry out appropriate confirmatory genetic tests.

Example of screening tests:

Case 1



Immunoglobulin profile

Total Immunoglobulin levels

| | | |
|---------------------------------------|-----------|---|
| Immunoglobulin-A RATE NEPHELOMETRY | <27 mg/dl | Neonate (4 Days) 0-2.2 Less than 2 yrs 14-108 Less than 12 yrs 29-270 Adult (20 -60 yrs) 70-400 >60 Yrs 90-410 |
| Immunoglobulin-G RATE NEPHELOMETRY | 361 mg/dL | New born - (4 days) 700-1480 Less than 2 yrs 500-1200 Less than 12 yrs 700-1650 Adult (20 -60 yrs) 700-1600 >60 Yrs - 600-1560 |
| Immunoglobulin-M RATE NEPHELOMETRY | 41 mg/dl | New born - (4 days) 5-30 Less than 2 yrs 43-239 Less than 12 yrs 50-260 Adult (20 -60 yrs) 40-230 >60 Yrs - 30-360 |

TBNK flow cytometry :


| Test Name TBNK (Flowcytometry) | Test Result | Biological Reference range |
|---|-------------------|--|
| TBNK- ABSOLUTE LYMPHOCYTE COUNT FLOWCYTOMETRY | 4887 cells/microL | 0 to 2 years : 700-11,9000 2 to 5 years : 1700 - 6900 5 to 10 years : 1100- 5900 10 to 16 years : 100 - 5300 Adults : 1000-2800 cells/microL |
| CD3 + : FLOWCYTOMETRY | 92.8% | % |
| CD19 + : FLOWCYTOMETRY | 1.0% | % |
| CD3 - (CD56 + CD16 +) FLOWCYTOMETRY | 5.2% | % |
| CD3 + (Absolute counts) FLOWCYTOMETRY | 4535 cells/microL | 0 to 2 years : 600-8000 2 to 5 years : 900 - 4500 5 to 10 years : 700- 4200 10 to 16 years : 800 - 3500 Adults : 700-2100 cells/microL |
| CD19 + (Absolute counts) FLOWCYTOMETRY | 49 cells/microL | 0 to 2 years : 600-8000 2 to 5 years : 900 - 4500 5 to 10 years : 700- 4200 10 to 16 years : 800 - 3500 Adults : 700-2100 cells/microL |
| CD3 - CD56 - CD16 + (Absolute counts) FLOWCYTOMETRY | 245 cells/microL | 0 to 2 years : 100-1400 2 to 5 years : 100 - 1000 5 to 10 years : 90- 900 10 to 16 years : 70 - 1200 Adults : 90-600 cells/microL |

Interpretation : "In view of decreased immunoglobulins and low B cell counts, Possibility of X linked agammaglobulinemia needs to considered

Case 2

- ▶ A 2-year old boy presented with repeated episodes of suppurative lymphadenitis. He was admitted with severe pneumonia.
- ▶ Blood investigations showed marked neutrophilic leukocytosis and thrombocytosis. He was suspected to have chronic granulomatous disease and evaluated further.

DHR test: A transport control sample was sent from an unrelated donor.




Legend: A and B are graphs from a transport control sample. A clear shift is seen in the neutrophil fluorescence post PMA stimulation, indicated by the blue arrow. In C and D which represent the test sample, there is hardly any shift in the fluorescence post PMA stimulation.

Interpretation: "Neutrophil oxidative burst is markedly reduced in comparison with the control sample. In view of history and clinical features, consistent with **chronic granulomatous disease**"

Case 3

- ▶ A 9-month-old girl with chronic diarrhea and pyoderma gangrenosum. She had a delayed cord fall at birth.
- ▶ Initial work-up showed marked neutrophilic leukocytosis. Clinical suspicion of Leucocyte Adhesion defect and Chronic granulomatous disease.

DHR test:



CD18 Flow cytometry for Leucocyte Adhesion defect:

CD18 expression on neutrophils - Control : 99.8%

CD18 expression on neutrophils - Test: 0.1%

Interpretation: "Almost complete absence of CD18 expression on neutrophils, consistent with Leucocyte Adhesion defect. Neutrophil respiratory burst appears preserved."

| Test/Profile | Components | Methodology | DoS Code | Sample Requirements | Reporting Time |
|--|---------------------------------------|----------------|----------|--|--|
| CD18 [For screening Leucocyte adhesion deficiency] | enumeration of CD18 | Flow cytometry | C0201 | 3 ml of whole blood in EDTA tube (Lavender top) & Sodium heparin. Ship refrigerated | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Dihdrorhdamine (DHR) | oxidative burst | Flow cytometry | D0029 | 3 ml (2 ml min.) whole blood in a 1 Green Top (Sodium Heparin) tube, ship refrigerated. DO NOT FREEZE | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Double Negative T- cells (DNT) | Analysis of double negative T - cells | Flow cytometry | D0031 | Collect 3 ml whole blood in EDTA tube (Lavender top). Ship refrigerated | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| CD 19 | CD 19 | Flow cytometry | F0088 | 3 mL (2 mL minimum) whole blood in 1 EDTA tube (Lavender Top) and 3 mL (2 mL minimum) whole blood in 1 Sodium Heparin tube (Green Top) or 2 mL (1 mL minimum) heparinized Bone marrow. Ship Immediately at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| CD 20 | CD 20 | Flow cytometry | F0091 | 3 mL (2 mL minimum) whole blood in 1 EDTA tube (Lavender Top) and 3 mL (2 mL minimum) whole blood in 1 Sodium Heparin tube (Green Top) or 2 mL (1 mL minimum) heparinized Bone marrow. Ship Immediately at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| CD4/CD8 absolute counts | CD3/4/8 & 45 | Flow cytometry | 10001 | Sample to be shipped at room temperature 48 hrs of collection | Sample reaching the lab before 2 PM on weekdays will be reported on the same day by 7 PM |

| Test/Profile | Components | Methodology | DoS Code | Sample Requirements | Reporting Time |
|--|---|-------------------|----------|---|--|
| Immuno-globulin IgA, Serum | IgA in Serum | Nephelometry | 10118 | 2 mL (1 mL min.) serum from 1 SST. Fresh sample preferred. Ship refrigerated. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Immuno-globulin IgE, Serum | IgE in Serum | Chemiluminescence | 10119 | 2 mL (1 mL min.) serum from 1 SST. Fresh sample preferred. Ship refrigerated. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Immuno-globulin IgG, Serum | IgG in Serum | Nephelometry | 10120 | 2 mL (1 mL min.) serum from 1 SST. Fresh sample preferred. Ship refrigerated. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Immuno-globulin IgM, Serum | IgM in Serum | Nephelometry | 10121 | 2 mL (1 mL min.) serum from 1 SST. Fresh sample preferred. Ship refrigerated. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Nitroblue Tetrazolium (NBT) | Tetrazolium nbt | Manual | NO022 | 3 mL (2 mL min.) Whole blood in 1 Green Top (Sodium Heparin) tube, ship refrigerated. DO NOT FREEZE | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| Primary immuno deficiency disorders (PIDs) Panel | NBT, TBNK, IgG, IgA, IgM, IgE, DHR, CD3, CD4, Cd8 | Flow cytometry | PO087 | For Patient- 2 mL (1 mL min.) serum from 1 SST & 2 mL (1 mL min) blood from 1 Green Top (Sodium Heparin) & 2 mL (1 mL min) Whole Blood Lavender Top (EDTA) tube. For control - 2ml whole blood in green top (Sodium heparin) from a normal control preferably unrelated. Ship refrigerated | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |
| TBNK | T Cell, B Cell & Natural Killer Cell enumeration | Flow cytometry | T0066 | 3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship immediately at 18 - 22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. | Sample reaching the lab before 1 PM on weekdays will be reported on the same day by 7 PM |

PARTNERS IN HEALTH



DR. SUJAY PRASAD

Hepatopathology,
Lymphoreticular system, Flow Cytometry
Technical Director
drsujay@neubergdiagnostics.com



DR. ANANTHVIKAS JAYARAM

Lymphoreticular, Hematooncology,
Flow Cytometry, GI Pathology,
Molecular Pathology, PID,
Hemoglobinopathies
Chief of Lab
ananthvikas@neuberganand.com



DR. PRADEEP KUMAR

PID, Hemoglobinopathies,
Preanalytics
pradeep_kv@neuberganand.com

FOR MORE DETAILS, CONTACT US AT



18004251974

neubergdiagnostics.com



Neuberg
DIAGNOSTICS

• India • UAE • South Africa • USA